



(Temperature self-adaption adjustment system)

KaniStar UC-2005/2007/2008 TSaaS ICU

TSaaS-ICU

ICU with Temperature self-adaption adjustment system

ICU- A tool used to improve quality of life for patients around the world. Kanistar is here to redefine the meaning of that care with our new interactive and innovative ICU line.


The purpose of any animal monitoring cage is to provide optimal care in hopes of aiding and speeding up recovery time for animal patients.

Kanistar ICU line does just that. At the machines most basic functions, it creates an environment of stable temperature and humidity, helping to provide a comfortable environment for the animal. These basic functions are great for newborn dogs and cats, especially born by caesarean section, as they cannot

regulate their own body temperature and are confronted by a high mortality rate caused by hypothermia. A warm and comfortably controlled environment helps to increase their survival rate. With advanced diagnosis, monitoring and treatment equipment and technology, our ICU helps to observe the qualitative and quantitative condition of your patient continuously and dynamically. The typical animal ICU machine on the market fails to close the gap between what we know human ICUs can do. In our redefining of the animal ICU, we hope to close the gap from being a simple enclosed space that provides temperature and humidity settings and offering a machine that can support multiple medical treatments, such as infusions and nebulization of medication, an oxygen therapy enclosure, the ability to purify the ambient air inside the containment, and more. Enjoy the animal ICU, redefined, in Kanistar's TSAAS intensive care units.



Features and Advantages



- 1. Infrared Physiotherapy
- 2. Built-in Medical Nebulization
- 3. Oxygen Concentration Control
- 4. Humidity Control (External Humidifier Needed)
- 5. Overhead Double Circulation Heating System
- 6. Adjustable Amber Light for Comfortable Illumination
- 7. White Examination Light
- 8. Self-Adaptive Temperature Control System (PRO Model Only)
- 9. IP Camera for Real-time Observation on Smart Devices (Pro Model&Optional)
- 10. Large Transparent Door to Easily Check Animal Status
- 11. Double Insulated Body for Better Heat and Oxygen Preservation
- 12. Triple Air Purification System Ultraviolet, Negative Ion, and External Circulation
- 13. Thermostatic Bed with Mutual Accurate Temperature Control
- 14. Aluminum frame on Front Door for Better Oxygen Sealing
- 15. Food/Medical Grade Stainless Steel Inner Body for Effective Sanitation Management
- 16. Large Touch Screen for Ease of Use and to Aide in More Visible Alerts
- 17. AI (Artificial Intelligence) Controlled Alarm System for Machine Errors like Abnormal Temp inside the Cabin, Abnormal Humidity, Abnormal Oxygen Concentration, and Abnormal Temp of the Thermostatic Bed
- 18. Hand Access for Easy Treatment with Ventilation Adjustment
- 20. Programmed Power Supply for External Devices like Humidifier
- 21. Uninterrupted Power Supply for Other Add-on Devices like Vital Sign Monitor
- 22. Oxygen Sealing Transfusion Portal
- 23. Infusion Post Fixture(optional)
- 24. Control Interface of Oxygen Generator
- 25. Easy to Stack Design, Two Levels for Better Space Efficiency.
- 26. Stack-on Detachable Wheel Base for Mobility
- 28. Interactive Indicator Light Band at the Front Bottom

Features and Advantages






- 19. Enclosing Large Enough to Accommodate a Large Dog Standing Comfortably
- 29. Power Switch of Air Conditioner
- 30. Port for Compressor of Nebulizer
- 31. Humidifier Port
- 32. Flowmeter of Oxygen
- 33. Oxygen Port
- 34. Air Conditioner

Dual drive constant temperature

Overhead double-duct indirect-blow circulating heating system

Indirect-blow refrigeration compressor air-conditioner





Physical cooling of the bottom deck (2023 optional)



Thermostatic bed with evenly-distrubuted heating and temperature control

Bidirectional control constant humidity



Breeding environment


Ultrasonic humidification technology Up to 75% high humidity at high temperature of 35 °C



Medical environment


Efficient dehumidification performance The humidity can be as low as 40%






High oxygen concentration control

- Intelligent stackable tower type oxygen concentrator with automatic control interface
- Control technology for balanced distribution in two/four compartments from single-side oxygen supply
- Intelligent Overlay Control Technology for Double/Four Cabins of two-side Oxygen Supply
- Conventional oxygen chamber oxygen supply (quickly reaching 50% oxygen concentration)
- High oxygen seal patch which helps to achieve up to 80% oxygen concentration





Complementary medicine

Built-in Nebulizer IR Infrared therapy technology

Available for atomization of the original solution without dilution; The atomized particles reach 0.5 microns and directly enter the lower respiratory tract for treatment.

Double oxygen supply holes/catheter holes/ IV pole



IR Infrared therapy technology

Enhance the vitality and regeneration ability of histiocyte, improve the metabolism of pets, diminish inflammation, sterilize, and accelerate wound healing







UV disinfection



The first: With bactericidal UVC ultraviolet light (wavelength 200~280 nm), Ultraviolet disinfection destroys the molecular structure of microbial cell DNA, to achieve the effect of sterilization.

High concentration negative ion



The second: High-concentration negative ion generators in a closed treatment chamber can significantly reduce odor pollution and enhance the immunity of animals.

Internal/external cycle ventilation control



The third: When oxygen supply is not required, the "external cycle" function can be on to improve the cabin air quality.

Double door and double cabin control



ICU Dual Cabin of fully Separated VersionDual cabin dual system of independent control



ICU Dual Cabin of Shared Version Shared control system for Dual Cabin



Isolation device-urine barrier for double cabin
Prevent the two compartments from being
soaked with pet urine




Prevent the two compartments from being soaked with pet urine

Prevent the two compartments from being soaked with pet urine

Lighting assistance

TSAAS ICU Creating a Rehabilitation Environment

Suitable temperature and humidity, continuous monitoring of body temperature, purification of sterile cabin, maternal quiet, soft illumination, exam lighting, high concentration oxygen chamber



White exam lamp Creating an examination environment for observing changes in diseases



Safety guarantee



Air conditioning intake to prevent hair suction

Easy to remove, built-in filter to avoid blockage or cross infection


Dual backup intelligent emergency vent

The emergency vent will open for the safety of pets' lives if there is oxygen sufficient or power failure in the cabin.

Air filter for cleaning

Top mounted air outlet filter with magnetic cover for easy replacement and cleaning


Medical care



Hand Access for Easy Treatment with Ventilation Adjustment

No need to open the cabin door, convenient for nursing operations Ventilation rotary cover can be rotated for rapid

ventilation of the cabin.







4-Level nursing level indicator

intensive care postoperative care hospitalized observation

mild care





Type 2007: Single cabin central touch screen



Type 2005: Dual cabin independent system central control touch screen



Type 2008: Dual cabin independent central control touch screen

Intelligent central control with large touch screen display

Simple and clear interface, with large digital display and touch control for various indicators, easy to identify and operate

TSAAS APP



Mobile phone Manage all ICUs and patients in the hospital (remote monitoring and management system)

Mobile phone Fully digital communication with pet owners anytime and anywhere on medical progress







Innovative IoT technology



Scan WeChat Open the TSAAS mini program



Remote management of ICU and sick pets

Remote monitoring and control of the status of the ICU and the temperature of pets, marking the ICU as a management device for the hospital. A doctor can efficiently manage all inpatient ICUs and pets in the hospital with just one mobile phone, and provide timely and accurate adjustment for pet conditions.

Continuous monitoring of pet temperature

7X24 hour Real time monitoring of changes in animal body temperature, alarming and automatic adjustment, enabling the ICU to iterate from a traditional incubator to a truly intelligent monitoring medical device.

