

42023 Surgical Cap, Miss

Surgical Cap

Product details

Size : One Size
Descriptive feature : Elastic, Extra comfort
Colour : Green
Sterile : Non-sterile
Brand : BARRIER

Images



Delivered items

42023-01

Sales released in: Australia, Austria, Belgium, Cyprus, Czechia, Denmark, Finland, Germany, Greece, Hungary, India, Ireland, Italy, Malaysia, Netherlands, New Zealand, Norway, Poland, Singapore, Slovakia, Slovenia, Sweden, Switzerland, United Kingdom of Great Britain and Northern Ireland

Shelf life: 5 years

Sterilization method: Non-sterile

Production Responsibility: WARMIŃSKA Sp. z o.o., ul. Wysokiej Bramy 31, 11-100 Lidzbark Warmiński, Poland

Packing information: First packaging layer is a cardboard dispenser box. Second layer is a corrugated board transport box.

Is suitable for Tray: No

Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Piece	1	7332430110480			
Dispenser box	120	7323190031078			

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.



Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Transport box	960	7323190031061	260x588x429 mm	65.6 dm ³	6.4 / 4.8 kg
Pallet	23040	7323190031054	800x1200x1850 mm		

Technical

Dimension

Dimension text	Dimension value
Fully stretched	290 mm
Unstretched	150 mm

Disposal instructions

Non-hazardous waste used BARRIER products and sterility barriers should, in the majority of cases, be classified as non-hazardous waste. They contain high amounts of energy and are well suited for incineration. BARRIER products do not contain any hazardous substances that can leach out if the products are land filled. Transport boxes are designed to fit existing recovery systems. The new BARRIER packaging system complies with the Packaging Waste Directive of the European Union.

Storage instructions

Mölnlycke Health Care recommends that BARRIER products are stored under normal storage conditions. All layers of packaging should be kept intact until access to the underlying layers is needed. Storage facilities for products only protected by the sterility barrier should be kept under conditions where low level of particulate air contamination prevail, so that it would not constitute a risk to the patient when the package is opened and the product is used.

Classifications

Regulation type(s)	MDR Class I ns	Locally Regulated	Locally Regulated	Locally Regulated
Intended Purpose :	The intended use of headwear is infection control by reducing the risk of transferring bacteria, hair or particles into the surgical wound	The intended use of headwear is infection control by reducing the risk of transferring bacteria, hair or particles into the surgical wound	The intended use of headwear is infection control by reducing the risk of transferring bacteria, hair or particles into the surgical wound	The intended use of headwear is infection control by reducing the risk of transferring bacteria, hair or particles into the surgical wound
MDR Classification Rule :	1			
Conformity Annexes :	IV			
Measuring Function :	No			
Notified body medical devices/PPE :	MPA	Local Authority New Zealand	Local Authority India	Local Authority Singapore

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.



Classifications

Regulation type(s)	Locally Regulated	AU Class I ns
Intended Purpose :	The intended use of headwear is infection control by reducing the risk of transferring bacteria, hair or particles into the surgical wound	The intended use of headwear is infection control by reducing the risk of transferring bacteria, hair or particles into the surgical wound
MDR Classification Rule :		
Conformity Annexes :		
Measuring Function :		
Notified body medical devices/PPE :	Local Authority Malaysia	TGA (Australia)

Applied standards : The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 62366, EN ISO 10993-5, EN ISO 10993-1, EN ISO 13485, EN ISO 9001, EN ISO 10993-10, ISO 15223-2, EN ISO 15223-1, ISO 14001

GMDN Code (Global Medical Device Nomenclature)

32297 Surgical cap, single-use

UNSPSC

42131611 Surgeon caps or hoods

Commodity Code

6505009090 Caps and hoods

CE Responsibility

Mölnlycke Health Care AB, Gamlestadsvägen 3c, Box 13080, SE-402 52 Goteborg Sweden

Basic UDIDI

73324300000000230JT

Find out more at www.molnlycke.com

