


Embryo transfer with 0.5 ml straws

Equine ET sheath , individually packaged, sterilized, 10/package front end opening [1]	19290/1050
side openings [2]	19290/1060
Stiff stylet , with white thumb rest	19290/1049
Flexible stylet , length: 65 cm	17209/1065
Finger support and ring , used with equine ET sheath and AI stylet [3]	19290/1048
Equine ET pipette , metal tip with side opening, flexible, disposable, sterilized, 1/package [4]	19290/1065
Stiff stylet , with white thumb rest	19290/1049
Flexible stylet , length: 65 cm	17209/1065
ET cannula , used with 0.5 ml straws, 1 side outlet, length: 55 cm, screw-off tip, stainless steel [5]	19290/0000



Embryo transfer with 0.25 ml straws

Deep intrauterine AI/ET set , introduction stylet, flexible AI/ET stylet, finger support, two fixation rings, for embryo deposition deep inside the uterine horn using 0.25 ml straws and DIUI sheaths	17015/0100
DIUI sheath for 0.25 ml straw	
10/bag	17015/0650
1/bag, sterilized	17015/0660



Embryo transfer > day 8

Large diameter pipette for equine AI and ET > day 8 , rounded ends, fits Luer syringe, length: 65 cm, sterilized	
1/bag	19290/1101
25/bag	19290/1105

