



M1001

THE GLOBAL NAME IN VETERINARY DENTISTRY

DENTAL MAT

Machine Washable at 60 °C / 140 °F

- Do not exceed 60°C / 140 °F
- Use standard detergent
- Recommend using a delicates washing net
- Low Spin (**High spin may cause damage to mat**)

Not suitable for:

- Machine drying
- Ironing
- Autoclaving
- High spin washing machine setting

Recycle at end of life – iM3 order code M1001

Made of Polyphenylene Ether (PPE) Plastic

Made in Germany for iM3.

iM3 Dental Ltd. Ireland - sales@im3dental.com
iM3 Pty Ltd. Australia - sales@im3vet.com
iM3 Inc. Vancouver, WA USA - info@im3usa.com



60°C/140°F

