Long Term Polyurethane Catheter Kits

INDICATIONS:

The Long Term MILACATH™ is indicated for venous administration of drugs or fluids, blood sampling, TPN and blood pressure monitoring. **Optional silver zirconium phosphate throughout the catheter body** including the inner and outer surface of the catheter is designed to inhibit bacteria growth without bacterial resistance.

RECOMMENDED INDWELL TIME:

1-30 days indwelling. The catheter is non-kinking and non-thrombogenic.

SELECT FROM THESE FEATURES TO DETERMINE ITEM NUMBER NEEDED ON PAGES 6-9:

- 1) Insertion Technique Guidewire Introduction or Peel Away Introduction.
- 2) Animal Size Large Animal/Equine or Small Animal/Dogs and Cats.
- 3) Number of Lumens and Length Single, Double, Triple or Quad Lumens at various lengths.
- 4) Catheter Type Standard or Antimicrobial.
- 5) Preferred Extensions Pre-Attached Extension Lines or Non-Attached Extension Lines.

All kits contain: Catheter, scalpel, introducer depending on insertion technique, needle free ports, suture wings, and gauze. Those that have extensions include clamps.

SMALL ANIMAL KITS include a fenestrated drape and 3-0 nylon suture.

LARGE ANIMAL KITS include a larger introducer and 2-0 nylon suture.

Single Lumen Catheters - Available with preattached or non-attached extension tube. The 14Ga catheter has a side and distal exit port.



Multi-Lumen Catheters - Separate lumens with separate exit points.









