

Test	Animal	Test Time	Sample Volume	Sample Type	Detection Range
<b>Cortisol</b>	Canine or Feline	15 min	75µl	Serum or Plasma (Lithium Heparin)	10-1000 nmol/L
<b>T4</b>	Canine or Feline	15 min	75µl	Serum or Plasma (Lithium Heparin)	8-100 nmo/L
<b>TSH</b>	Canine or Feline	15 min	75µl	Serum or Plasma (Lithium Heparin)	0,1-100 U/L
<b>Progesterone</b>	Canine	15 min	75µl	Serum or Plasma (Lithium Heparin)	1-50 ng/ml
<b>cCRP</b>	Canine	3 min	10µl	Serum or Plasma (Lithium Heparin)	2-250 mg/L
<b>fSAA</b>	Feline	3 min	10µl	Serum or Plasma (Lithium Heparin)	5-300 mg/L
<b>D-Dimer</b>	Canine or Feline	5 min	10µl	Plasma (Lithium Heparin)	50-10,000 ng/ml
<b>HbA1c</b>	Canine	5 min	10µl	Whole Blood (EDTA)	4-14,5%
<b>cPL</b>	Canine	15 min	75µl	Serum or Plasma (Lithium Heparin)	10-1000 µg/L
<b>fPL</b>	Feline	15 min	75 µl	Serum or Plasma (Lithium Heparin)	0,5-100 µg/L
<b>NT-proBNP</b>	Canine	15 min	75 µl	Serum or Plasma (Lithium Heparin)	5-5000 pmol/L