

NEW PRODUCT



IV Break-Away Device

Connect to IV administration line and prevent IV catheter from being dislodged



- No more T-ports snapped or IVs dislodged!
- Lightweight and low-profile: QUICK and EASY to add to IV lines
- No blood loss or mess from broken IV lines
- No more wasted time, money, and stress associated with restarting IVs in unruly patients!
- Supports gravity flow rates of >220mL/min
- Easy to reassemble



How it works

Simply
disinfect
and
reconnect!

Each end AUTOMATICALLY
closes when disconnected!



Item #IVBA

