

Safety-Lancet Product Range

					
Version	Mini	Normal	Extra	Super	Neonatal
Order no.	85.1015	85.1016	85.1017	85.1018	85.1019
Needle size	28 G	21 G	18 G	Blade 1.5 mm	Blade 1.5 mm
Penetration depth	1.6 mm	1.8 mm	1.8 mm	1.6 mm	1.2 mm
Blood volume app.* after 20 sec after 120 sec	50 µl 150 µl	100 µl 500 µl	150 µl 750 µl	150 µl 900 µl	100 µl 450 µl
Puncture site	finger, ear, heel	finger, ear, heel	finger, ear, heel	finger, ear, heel	finger, heel

* According to study with 30 adult test persons.

Please be aware that the selection of the appropriate lancet depends upon various factors, such as skin type, blood flow to the puncture site and the correct application of the Safety-Lancet onto the skin. In case of doubt, it is recommended to test a selection of different Safety-Lancets prior to using the products for routine applications. The information given above is a recommendation only and does not, under any circumstances replace medical, scientific or technological advice.