

<b>REF 19E-370.01</b>	<b>Product</b> container filter with indicator	<b>Size</b> 190 mm
<b>Product structure:</b>	sterilization crepe paper white, 60 g/m <sup>2</sup> according to DIN EN ISO 11607-1	
<b>Performance:</b>	round with diameter = 190 mm, with tab	
<b>Indicator:</b>	steam indicator Type 1 according to DIN EN ISO 11140-1:2014 colour change from pink to dark brown - black	
<b>Printing colour:</b>	1 – colour, black, sterilizable	
<b>Toxicity:</b>	toxic effect is not known	
<b>Product description:</b>	Container filters made of paper for single use for Aesculap-Containers Round filters with steam indicator and a perforated hole The Container filters are for single use and not recyclable.	
<b>Instruction for use:</b>	Place one filter over each perforated section on the inside of the SterilContainer lid. Secure each filter with the filter retention plate.	

### Compatibility

#### reusable sterilization containers:

Aesculap:	<input checked="" type="checkbox"/> suitable	<input type="checkbox"/> not suitable
Wagner:	<input type="checkbox"/> suitable	<input checked="" type="checkbox"/> not suitable
Martin:	<input type="checkbox"/> suitable	<input checked="" type="checkbox"/> not suitable
Medicon:	<input type="checkbox"/> suitable	<input checked="" type="checkbox"/> not suitable

### Sterilisation prozess:

STEAM	<input checked="" type="checkbox"/> suitable	<input type="checkbox"/> not suitable
EO	<input type="checkbox"/> suitable	<input checked="" type="checkbox"/> not suitable
FORM	<input type="checkbox"/> suitable	<input checked="" type="checkbox"/> not suitable
VH2O2	<input type="checkbox"/> suitable	<input checked="" type="checkbox"/> not suitable

### Layout:

