

## Standard and Primary Sets

### VL ST00

ref. M46441000S

#### Approximate set length:

285 cm


#### Approximate priming volume:

25 ml

#### Main materials:

PVC (tubing: diameter 3.0 x 4.1 mm)

Silicone (tubing: diameter 3.41 x 5.11 mm)



#### Advantages

- Anti free-flow safety clamp
- DEHP as plasticizer is not used
- Latex-free

#### Packaging

- Case quantity: 30 x 1

#### Typical applications

- Infusion of drugs / fluids
  - having no particular filtration requirements
  - available in bags or bottles
  - having no particular incompatibility

#### Not to be used for / in

- Blood cells
- Infusion of drugs known to be incompatible with PVC (e.g trinitrine, isosorbide dinitrate, sodium nitroprussiate...)
- Light sensitive drugs (unless other precautions are taken to protect the drug against direct light)
- Primary / secondary mode with Volumat MC Agilia / Agilia VP MC

Due to our policy of continuous product improvement, as well as changes in standards, the features described in this document may change. Please consult us for the most current information.