

Ear Diagnostic Set

Diagnostic otoscopy in veterinary medicine is one of the most important fields of application in endoscopy. Disorders of the external ear are common in dogs, therefore, the easy accessibility afforded by this technique makes otoendoscopy an ideal examination method.

A rigid endoscope can also be utilized to examine the middle ear.

KARL STORZ now offers a new veterinary otoscopy set which consists of an otoscope with a 0° forward-oblique telescope and a 3-way stopcock.

The shape of the telescope and its length facilitate thorough examination even in large breed dogs.

The instrument set can also be used for treatment, i.e. flushing the ear canal, suction, removal of foreign bodies, etc.

Used in conjunction with a KARL STORZ camera system and a video recorder or printer, it provides a precise means of assessing clinical results and documenting findings.



67260 OSA

67260 OSA **Veterinary Otoscope Set** including:

Otoscope, for dogs and cats, **autoclavable**, diameter 5 mm, length 8.5 cm, with working channel diameter 5 Fr.

Stopcock Attachment, with integrated working channel and LUER-Lock adaptor for irrigation

Optional Accessories



60651 BK

Cleaning Brush, flexible, for use with working channels up to 2.3 mm, length 25 cm



NEW 69003 HH

Otoscope Camera Head Holder, for combining Otoscope 67260 AA with **non-autoclavable** TELECAM camera head and Veterinary Video Camera III camera head



NEW 15006 D

"ULTRA STOP" Antifog Solution, 15 ml, atomizer bottle

Optional accessories for ear flushing and suction see chapter 11, UNITS, VETPUMP® II

Camera systems and accessories for documentation see chapter 12, EXCERPT FROM TELEPRESENCE

Components/Spare Parts see chapter 14